



Product Specification

Product name	Lolina® Osteoblast Mineralization Supplement
Cat.No.	NaC201201
Size	5ml
Storage and shipping	Store ObMS at -20°C before adding to basal medium. Dry ice.

Product Description

Lolina® Osteoblast Mineralization Supplement (ObMS) is a medium supplement designed for mineralization of primary osteoblasts as determined by Alizarin Red staining in vitro. It is a sterile, concentrated (100X) solution which contains growth factors, hormones, and proteins necessary for osteoblast differentiation. The supplement is designed as an additive for Osteoblast Mineralization Medium (ObMM) and should be used in conjunction with that medium.

Prepare for use

Thaw ObMS at 37°C. Gently tilt the tube several times during thawing to help the contents dissolve. Make sure the contents of the supplement are completely dissolved into solution before adding to the medium. Rinse the bottle and tubes with 70% ethanol, and then wipe to remove excess. Remove the cap, being careful not to touch the interior threads with fingers. Add ObMS and other components (FBS and P/S solution) into the basal medium in a sterile field. Mix well and the reconstituted medium is ready for use. Since several components of ObMM are light-labile, it is recommended that the medium not be exposed to light for lengthy periods of time. If the medium is warmed prior to use, do not exceed 37°C. When stored in the dark at 4°C, the reconstituted medium is stable for one month.

Caution: If handled improperly, some components of the medium may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear. Dispose of properly.

Note

ObMS is for research use only. It is not approved for human or animal use, or for application in in vitro diagnostic procedures.