



Product Specification

Product name	Lolina® Trabecular Meshwork Cell Growth Supplement
Cat.No.	NaC127101
Size	5ml
Storage and shipping	Store TMC GS at -20°C; once added to medium, store at 4°C and protect from light. Do not refreeze after thawing. Dry ice.

Product Description

Lolina® Trabecular Meshwork Cell Growth Supplement (TMC GS) is a sterile, concentrated (100X) solution which contains growth factors, hormones, and proteins necessary for the culture of normal trabecular meshwork cells. The supplement is formulated (quantitatively and qualitatively) to provide an optimally balanced environment that selectively supports the growth of normal trabecular meshwork cells in vitro.

Instructions for Use

For best results, use TMC GS in conjunction with 500 ml Trabecular Meshwork Cell Medium-b (TMC M-b) and 10 ml fetal bovine serum (FBS).

Thaw TMC GS at 37°C. Gently tilt the tube several times to ensure complete mixing. Spray the medium bottle and supplement tube with 70% ethanol and wipe to remove excess liquid. In a sterile field, remove the caps without touching the interior threads with fingers. Add TMC GS to the basal medium and mix well. Since several components are light-labile, the medium should not be exposed to light for extended periods. We do not recommend warming medium in a 37°C water bath prior to use. When stored in the dark at 4°C, the reconstituted medium is stable for one month.

Caution: If handled improperly, some components of the medium may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear. Dispose of properly.

Note

TMC GS is for research use only. It is not approved for human or animal use, or for application in in vitro diagnostic procedures.